

# POCKET SHARP CONTAINER



|               |                   |
|---------------|-------------------|
| Product Code  | 25780             |
| Unit of sale  | box of 70 pcs     |
| Minimum order | 1                 |
| Type          | No medical device |
| Class         |                   |

|         |               |
|---------|---------------|
| UK-REP  | Not required  |
| CH-REP  | Not required  |
| EAN/UPC | 8013183134143 |
| PARAF   | 977810833     |

## Description

POCKET SHARP CONTAINER 0.4 l - 12.3x6.2xh 16.5 cm

Single use container for both sharp and cutting waste disposal.

The permanent and inviolable closure that prevents waste leakage and the high container resistance guarantee top quality and safety.

Equipped with a patented double opening system, it can be used by just one hand.

The double opening includes one side with needles and insulin pen removal, and the other side with a large opening for entire syringe disposal.

- Features:
- Pocket size
  - Ergonomic shape
  - Practical opening system
  - Needles removal
  - Insulin pen removal device
  - Maximum filling level
  - Shock resistant
  - Puncture resistant
  - Sturdy handle
  - Possible side fastening
  - Autoclavable at 134°C for 18'
  - No emission of any harmful gas during the incineration process

Made in Italy.