

DISPO SHARP CONTAINER - 0.5 l



Product Code	25781
Unit of sale	box of 187 pcs
Minimum order	1
Type	No medical device
Class	

UK-REP	Not required
CH-REP	Not required
EAN/UPC	8013183111786

Description

DISPO SHARP CONTAINER 0.5 l - 8.3x h 15.1 cm

Single use conical shape puncture resistant container for the disposal of sharp and cutting waste.

The permanent and inviolable closure that prevents waste leakage and the high container resistance guarantee top quality and safety.

Temporary and permanent closures with a wide opening for the disposal of larger materials, and three indents for the needles and insulin pen removal.

- Features:
- Pocket size
 - Practical opening system
 - Needles removal
 - Insulin pen removal device
 - Maximum filling level
 - Shock resistant
 - Puncture resistant
 - Sturdy handle
 - Possible side fastening
 - Autoclavable at 134°C for 18'
 - No emission of any harmful gas during the incineration process

Made in Italy.