

CS PLUS LINE SHARP CONTAINER - 10 I



Product Code	25787
Unit of sale	box of 33 pcs
Minimum order	1
Type	No medical device
Class	

UK-REP	Not required
CH-REP	Not required
EAN/UPC	8013183143596

Description

CS PLUS line sharps containers
• CS Plus 10 I - 22x22xh 32.5 cm

Full range of Sharps containers for the disposal of sharp or cutting waste.
Equipped with temporary and definitive closure, certified hermetic seal, do not release any harmful emission during incineration.
The containers hold the UN mark for the transportation of dangerous materials.

- practical opening system
- temporary and permanent closure
- needles removal
- insulin pen removal device
- maximum filling level
- shock resistant
- puncture resistant
- easy to use
- easy to assemble
- excellent stability
- square or round shaped
- sturdy handle for easier transport
- customisable label for quantities
- autoclavable at 134°C for 18'
- no emission of any harmful gas during the incineration process

Made in Italy.