

PBS LINE SHARP CONTAINER 12 I



Product Code	25795
Unit of sale	box of 25 pcs
Minimum order	1
Type	No medical device
Class	

UK-REP	Not required
CH-REP	Not required

EAN/UPC	8013183111991
PARAF	977810910

Description

PBS sharps containers

- PBS 12 I - Ø 25xh 31.3 cm

Full range of Sharps containers for the disposal of sharp or cutting waste.
Equipped with temporary and definitive closure, certified hermetic seal, do not release any harmful emission during incineration.
The containers hold the UN mark for the transportation of dangerous materials.

- practical opening system
- temporary and permanent closure
- needles removal
- insulin pen removal device
- maximum filling level
- shock resistant
- puncture resistant
- easy to use
- easy to assemble
- excellent stability
- square or round shaped
- sturdy handle for easier transport
- customisable label for quantities
- autoclavable at 134°C for 18'
- no emission of any harmful gas during the incineration process

Made in Italy.