

MAGNUM LONG SHARP CONTAINER 6.3 I



Product Code	25797
Unit of sale	1 pc
Minimum order	1
Type	No medical device
Class	

UK-REP	Not required
CH-REP	Not required

EAN/UPC	8013183134167
PARAF	977810922

Description

LONG SHARP CONTAINERS
• Magnum L 6.3 I - 17.2x12xh 61.7 cm

Disposable containers for the elimination of long objects up to 72 cm, such as biopsy forceps and sharps.
Large opening, for elimination of bulky objects, temporary and definitive closure, ergonomic handle.
Autoclavable at 134°C for 18'

Made in Italy.