

SINGLE PATIENT USE F.O. LED BATTERY HANDLE "RTU" (ready to use) 145 kLux - pediatric



Product Code	58018
Unit of sale	box of 10 pcs
Minimum order	1
Type	Medical device
Class	I S

RDM (NSIS)	2365941
CND	Z12021080
EAN/UPC	8961102788703

Description

SINGLE PATIENT F.O. LED BATTERY HANDLES 2.5 V

Battery handle "RTU" (ready to use) - pediatric - box of 10

- lightweight, suitable for everyday hospital use
- single patient use, prevent cross infection
- compatible with ISO 7376 Green System
- perfect solution to cleaning problems confronting EMS field
- Ready To Use (RTU) handles supplied with disposable eco-friendly preloaded batteries
- suitable also for reusable blades

Multilingual manual: GB, FR, IT, ES, PT, DE, RO, GR, Arabic.

Technical Specifications

CRI: >110
Material: ABS
Lux: 145K
Temperature: 4500K
Batteries: mercury free