


ELINE FLAT REELS 200 m x 50 mm AUTOCLAVE



Product Code	35803
Unit of sale	box of 8 pcs
Minimum order	1
Type	Medical device
Class	I

RDM (NSIS)	2755779
CND	S01020101
EAN/UPC	8023279358032
GMDN	13735

Description

FLAT REELS for sterilization by steam or ethylene oxide gas

- FLAT REELS
- External: 200 m x 50 mm
- Internal: 200 m x 36 mm

Made of transparent, polyester-polypropylene film coupled to medical grade paper which permit steam and gas penetration during sterilization.
Steam and gas indicators provide clear colour changes when exposed to sterilization process.

Made in Italy.