

LANCART TROLLEY 60x40xh80 cm



Product Code	45810
Unit of sale	1 pc
Minimum order	1
Type	No medical device
Class	

UK-REP	Not required
CH-REP	Not required
EAN/UPC	8023279458107
PARAF	951840622

Description

Lancart trolley - small

- Size: 60x40xh 80
- Shelves: 2
- Drawer: -

Stainless steel trolleys mounted on castors Ø 80 mm.
Maximum trolley load: 60 kg.
Range of optional accessories to customize your trolley.
Delivered in kit form.

Made in Italy.