



Product Code	45811
Unit of sale	1 pc
Minimum order	1
Type	No medical device
Class	

UK-REP	Not required
CH-REP	Not required
EAN/UPC	8023279458114

## Description

Lancart trolley - medium  
- Size: 72x50xh 80  
- Shelves: 2  
- Drawer: -

Stainless steel trolleys mounted on castors Ø 80 mm.  
Maximum trolley load: 60 kg.  
Range of optional accessories to customize your trolley.  
Delivered in kit form.

Made in Italy.