

LANCART TER TROLLEY 70x50xh80 cm - 3 shelves



Product Code 45813

Unit of sale 1 pc

Minimum order 1

Type No medical device

Class

UK-REP Not required

CH-REP Not required

EAN/UPC 8023279458138

Description

Lancart Ter trolley - medium
- Size: 72x50xh 80
- Shelves: 3
- Drawer: -

Stainless steel trolleys mounted on castors Ø 80 mm.
Maximum trolley load: 60 kg.
Range of optional accessories to customize your trolley.
Delivered in kit form.

Made in Italy.