

ELINE FLAT REELS 200 m x 400 mm AUTOCLAVE



Product Code	35816
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	Yes
CH-REP	Yes

RDM (NSIS)	2755779
CND	S01020101
EAN/UPC	8023279358162
PARAF	903118533
GMDN	13735

Description

FLAT REELS for sterilization by steam or ethylene oxide gas

- FLAT REELS
- External: 200 m x 400 mm
- Internal: 200 m x 378 mm

Made of transparent, polyester-polypropylene film coupled to medical grade paper which permit steam and gas penetration during sterilization.
Steam and gas indicators provide clear colour changes when exposed to sterilization process.

Made in Italy.