

DRY HEAT TAPE 19 mm - hot air sterilizer



Product Code	35842
Unit of sale	box of 12 pcs
Minimum order	1
Type	No medical device
Class	

UK-REP	Not required
CH-REP	Not required
EAN/UPC	8001814308575

Description

DRY HEAT TAPE INDICATORS

Tape indicator - 50 m: Dry Heat Tape - hot air sterilizer
Application: High temperature
Width 19 mm

Gima range of tapes is specially formulated for packing purposes as well as to show visibly that the pack has undergone a sterilization cycle.
The backing of the different tapes is made of strong, semicreped, impregnated paper which gives excellent conformability during the application.
Sterile tapes unwind easily and can be torn by hand without difficulty.
Special adhesives formulated to withstand high temperature and humidity ensure excellent bond to paper, plastic, non-wovens, board, metal, glass and linen.
Each tape is individually packed with due date.

Dry tape:
for heat cycles at 160° - minimum 40 minutes from green to brown.

Multilingual instructions: GB, FR, IT, ES, PT, DE, NL, DK, FI, RO, HR, GR.

Made in Italy.