



STEAM STAR INDICATOR - autoclave



Product Code	35849
Unit of sale	box of 1000 pcs
Minimum order	1
Type	No medical device
Class	
UK-REP	Not required
CH-REP	Not required
EAN/UPC	3760243050022

Description

Indicator for sterilization by autoclave with pressure 1 to 2 bar 120°C/135°C.
No special stocking requirements.
Safety and reliability.
Made in Japan.

Notice "Stérilisé" and three stars drawing change colour from pink to brown to show effectiveness of 121°C - 15 minutes or 134°C - 3 minutes and 30" sterilization and to black to show safety margin of 121°C - 20 minutes and 134°C - 5 minutes sterilization.