



Product Code	22585
Unit of sale	box of 10 pcs
Minimum order	1
Type	Medical device
Class	II A

UK-REP	Yes
CH-REP	No

RDM (NSIS)	1774717
CND	A010299
EAN/UPC	8023279225853
UMDNS	13-230
PARAF	977810389

Description

Sterile, ready-to-use, disposable punches eliminate need for sterilization.

- Biopsy Punches: Ø 3.5 mm
- Blade length: 8 mm

Seamless, stainless steel sharp cutting edge, ribbed handle for sure grip and control.

Size embossed on each punch for easy identification. Wide range of sizes.

Each punch is individually sealed in a peel-apart see-through package.

Clinical areas: dermatology, gynaecology, podiatry.

Multilanguage box and instructions: GB, FR, IT, ES, PT, DE, PL, GR, RU, FI, SE, NL, RO, BG, HR, HU, NO, Arab.

Seamless Stainless Steel Cutting Edge ensures a consistently smooth cutting surface, allowing optimal specimen quality.