

ELINE SELF-SEAL POUCHES 90x165 mm



Product Code	35857
Unit of sale	box of 2400 pcs
Minimum order	1
Type	Medical device
Class	I

UK-REP	Yes
CH-REP	Yes

RDM (NSIS)	2755769
CND	S01010101
EAN/UPC	8023279358575
GMDN	13735

Description

SELF-SEAL POUCHES FOR STERILIZATION BY STEAM OR ETO

- Self-seal (packed in box of 200 pouches inside cartoon)
- External (m x mm): 90 x 165
- Internal (m x mm): 68 x 150

Light green, polyester-polypropylene transparent film coupled to medical grade paper which permit steam and gas penetration during sterilization. Steam and gas indicators change colour when exposed to sterilization.

1. Reinforced corners to prevent dust from collecting near to the opening, reducing the risk of contamination.
2. The laminated transparent film complies with the German standard EN 10/204. The medical paper conforms to: DIN 58953 (German), Pharmacopoeia 589 (French), pr EN 868-3 (European) and BS 6256 (English).
3. The chevrons allow rapid and correct opening.
4. These consist of one layer of laminated transparent film (special type polyester and copolymer polypropylene) and white damp-resistant medical Kraft paper, specially designed to assure a high anti-bacterial protection.
5. The heat-sealed lines (Multiline) are highly resistant and spaced equally; the external welding runs along the entire length of the pouch.

Reels and pouches are certified and manufactured according to DIN 58953 BS 6871, EN 868-5, CN 13485.

Made in Italy.