

# OTOCLEAR EAR SPRAY WASH KIT - bottle 500 ml



|               |                |
|---------------|----------------|
| Product Code  | 25858          |
| Unit of sale  | 1 kit          |
| Minimum order | 1              |
| Type          | Medical device |
| Class         | I              |

|            |               |
|------------|---------------|
| RDM (NSIS) | 2681602       |
| CND        | Z12149099     |
| EAN/UPC    | 8023279258585 |

## Description

Otoclear ear spray wash kit

- manual, intermittent lavage
- delivers greater volume compared to syringes or other manual devices
- efficient design allows for two handed procedure

Improved patient comfort thanks to the soft "gentle touch" design  
Thermometer for water temperature

Kit includes:

- 1 spray wash bottle
- 20 otoclear tips
- 1 ear basin.

Made in the USA.