

OTOCLEAR AQUABOT EAR SPRAY WASH KIT - bottle 650 ml



Product Code	25859
Unit of sale	1 kit
Minimum order	1
Type	Medical device
Class	I

RDM (NSIS)	2681602
CND	Z12149099
GTIN	00859911004894

Description

OTOCLEAR AQUABOT® EAR SPRAY WASH - 650 ml

- provides continuous stream without manually pumping
- delivers greater volume compared to syringes or other manual devices
- efficient design allows for two-handed procedure

Kit includes: 1 spray wash bottle, 5 otoclear tips, 1 ear basin

Made in the USA.