

PILL DISPENSER with Buettooth



Product Code	25871
Unit of sale	1 pc
Minimum order	1
Type	No medical device
Class	

EAN/UPC	8023279258714
PARAF	950012385

Description

WIRELESS PILL DISPENSER

Automatic wireless pill dispenser to manage and record pills' dosage data with Healbox App and cloud server.

It lets patients and caregivers to program dosages up to a maximum of 28 consecutive days.

It also allows to give right dose to right patients at right time.

The device has audio and visual alerts up to 9 alarms per day, which sound for 30 minutes or until the dispenser is flipped.

It automatically rotates when it's time for the next dose.

28 compartments with 24 aspirin sized pills load capacity.

Lid can be locked with metal key.

For hospital, community service or personal use.

Supplied with 1 medication tray, 1 set of dosage rings, 1 metal key, 4xAA batteries and multilingual instructions (GB, FR, IT, ES, PT).

Easy setup only three buttons to setup time/alarm/volume

Bubble ring with cover and lock system for moisture-proof

Locking key included for secure lid

Technical Specifications

Compartments: 28 slots
Capacity: 24 Aspirin-sized pills
Material: ABS
LED display size: 60x30 mm
Bluetooth: 10 meters limited range
Battery: 4xAA batteries
Size: 225x h 68 mm
Weight: 600 g