



Product Code	58841
Unit of sale	box of 50 pcs
Minimum order	1
Type	Medical device
Class	I S

RDM (NSIS)	2514613
CND	M04010202
EAN/UPC	8015143008754
PARAF	987475340

Description

HYPOR FILM DRESSING FOR FIXATION - STERILE

- Hypor film 10 x 12.5 cm

Sterile adhesive dressing on transparent polyurethane film apt to fix catheters, gauzes, bandages and other dressings. The films acts as water and bacterial barrier.
Permeable to air and water vapor, lets the skin breathe and avoids the risk of maceration. Individually packed.
Extremely elastic substrate and highly conformable to sensitive areas.
Reduces the risk of fluid leaking and infection.

Multilingual box: GB, FR, IT, ES, DE.
Medical Device Class I.

Rounded edges to prevent the dressing from peeling off.
Hypoallergenic acrylic adhesive, dermatologically tested.
Protective liners easy removable with gloves.
Made in Italy.