


# ELINE FLAT POUCHES 250x400 mm



|               |                |
|---------------|----------------|
| Product Code  | 35888          |
| Unit of sale  | box of 500 pcs |
| Minimum order | 1              |
| Type          | Medical device |
| Class         | I              |

|        |     |
|--------|-----|
| UK-REP | Yes |
| CH-REP | Yes |

|            |               |
|------------|---------------|
| RDM (NSIS) | 2755772       |
| CND        | S01010101     |
| EAN/UPC    | 8023279358889 |
| GMDN       | 13735         |

## Description

### FLAT POUCHES FOR STERILIZATION BY STEAM OR ETO

- External (m x mm): 250 x 400
- Internal (m x mm): 228 x 385

Light green, polyester-polypropylene transparent film coupled to medical grade paper which permit steam and gas penetration during sterilization. Steam and gas indicators change colour when exposed to sterilization.

1. Reinforced corners to prevent dust from collecting near to the opening, reducing the risk of contamination.
2. The laminated transparent film complies with the German standard EN 10/204. The medical paper conforms to: DIN 58953 (German), Pharmacopoeia 589 (French), pr EN 868-3 (European) and BS 6256 (English).
3. The chevrons allow rapid and correct opening.
4. These consist of one layer of laminated transparent film (special type polyester and copolymer polypropylene) and white damp-resistant medical Kraft paper, specially designed to assure a high anti-bacterial protection.
5. The heat-sealed lines (Multiline) are highly resistant and spaced equally; the external welding runs along the entire length of the pouch.

Reels and pouches are certified and manufactured according to DIN 58953 BS 6871, EN 868-5, CN 13485.

Made in Italy.