

ELASTIC TUBULAR NETTING C/D for elbow, arm, foot - latex free



Product Code	58904
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

RDM (NSIS)	2272650
CND	M03030103
EAN/UPC	8033971550429

Description

Surgifix® elasticated tubular netting

- Material: Zerolatex (polyamide 86%, elastomer 14%)
- Range: C/D
- Net for elbow, arm, foot
- Length under tension: 25 m
- Width: 48 mm (Different width between Poly and Zerolatex is due to different elasticity of materials)

Surgifix® is a unique elastic tubular netting, light, soft and breathable.

It is the ideal solution for gauze or dressing retention.

Its unique structure makes it suitable for any part of the body, thus assuring perfect dressing adhesion under ideal compression grade and offering unrivalled patient's comfort.

Surgifix can be cut at the requested length.

Two ranges: Poly (polyester 38%, polyamide 35%, latex 27%) and Zerolatex (polyamide 86%, elastomer 14%).

Multilingual box (GB, FR, IT, ES, DE) and instructions (GB, FR, IT, ES, PT, DE, HU, PL, BG, GR, RO, SE, SI, NO, NL, LT, LV, CZ, DK, FI, RU, TR, SK, EE, JP, Arabic).

Made in Italy.