

BAG HOLDER TROLLEY foot operated - 3 bags



Product Code	45906
Unit of sale	1 pc
Minimum order	1
Type	No medical device
Class	

UK-REP	Not required
CH-REP	Not required
EAN/UPC	8023279459067

Description

BAG HOLDER TROLLEY - 3 bags

Bag holder trolleys with pedal operated opening.

Rigid, sturdy, self-supporting ABS plastic base with edges.

The bag holding ring (35x25.5 cm) is made of glass reinforced nylon with a rubber gasket to hold the bag in place.

Two vertical tubes made of AISI 304 stainless steel are welded to the base and carry the bag holding ring.

Non marking grey rubber swivelling wheels Ø 80 mm (two with pedal operated brake).

Pedal system equipped with shock absorber for silent lid closure.

Made in Italy.