


GIMA D-400 DIGITAL SEALING MACHINE without printer

230V



Product Code	35909
Unit of sale	1 pc
Minimum order	1
Type	No medical device
Class	

UK-REP	Not required
CH-REP	Not required

EAN/UPC	8023279359091
---------	---------------

Description

Continuous sealing machine for medical pouches according to DIN rules 58953-7:2010 without printer. It allow a 12 mm multi-line seal and features stand-by function and automatic infeed through photocell.

D-400 is a sophisticated unit controlled by a microprocessor and equipped with a device to adjust temperature (displayed by LCD) and printing parameters. The sealing pressure is automatically adjusted according to the material thickness. It also guarantees 30 mm safety distance between seal and product, according to DIN rules. Easy to use and maintenance-free, they are the ideal solution for clinics and hospitals. Compatible for sealing with Tyvek rolls and envelopes.

Made in Italy.

Technical Specifications

- Infeed speed: 6 m/min
- Temperature control: by digital thermoregulator
- Sealing temperature: 40° - 200° C
- Temperature tolerance: ±5%
- Seal: Multiline
- Seal width: 12 mm
- Adjustable free edge: 0 ÷ 30
- Power supply: 230V-50 Hz - single phase
- Absorption: 500 W
- Dimensions: 473 x 235 x 181 mm
- Net weight: 11,5 kg