



Product Code	59099
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
--------	----

CND	Z12139005
EAN/UPC	4058368096643
GMDN	16340

## Description

### HEBU GOLD II PLASTER SAW

The new oscillating professional saw can quickly, safely, quietly and accurately remove all plaster of Paris and synthetic dressings. High-performance motor with electronic speed regulation and a gradual start-up, gently on materials. Motor is cooled using the circulation principle. Cool air is suctioned in at the end of the motor, to protect power aggregate from resin dust and residues. Made of high quality materials.

Made in Germany.

## Technical Specifications

Voltage: 220-240 V, 50-60 Hz  
 Power: 250 W  
 Strokes: 12,000 - 21,000 oscillations/min  
 Length without cord: 340 mm  
 Weight with/without cable: 1.3/1.7 Kg  
 Noise: 68-82 dB  
 Cord length: 4.5 m  
 IP code: IP20

## Standard accessories

Saw blade Ø 50 mm  
 Saw blade Ø 65 mm  
 One pair of wrenches  
 User manual: GB, FR, IT, ES, DE