

LAUNDRY TROLLEY 3 shelves



Product Code	45913
Unit of sale	1 pc
Minimum order	1
Type	No medical device
Class	

UK-REP	No
CH-REP	No
EAN/UPC	8023279459135

Description

LAUNDRY TROLLEY

Laundry trolley with stainless steel structure.
Canvas bag and three shelves.
Silent and antistatic castors Ø 7.5 cm with brake assembled on the vertical column line for better stability.

Size 90x45xh 105 cm.
Weight: 14 kg.
Packing size: 93x50x97 cm
Supplied with canvas bag.

Instructions: GB, FR, IT, ES, PT, DE, NL, PL, RO, CZ, HU, SE, FI, HR, GR, BG, LT.

Made in Italy.