



Product Code	59211
Unit of sale	box of 10 pcs
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	No

RDM (NSIS)	2294827
CND	M0303010201
EAN/UPC	8023279592115

## Description

### ELASTOMULL® COHESIVE BANDAGES

Elastic retention bandage. It provides wound care splint retention and light support.  
Comfortable non woven fabric (cotton, viscose and polyamide) coated with synthetic latex.  
Stretched length 4 m.

- BSN code: 45471-00
- Elastomull® Haft 4 m x 6 cm
- Box of 10 pcs.
- Multilanguage box: GB, FR, IT, ES, DE, PT, NL, SE.