

ELASTOMULL HAFT COHESIVE BANDAGE 4 m x 10 cm



Product Code	59213
Unit of sale	box of 10 pcs
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	No

RDM (NSIS)	2294822
CND	M0303010201
EAN/UPC	4042809021912
PARAF	900161757

Description

ELASTOMULL® COHESIVE BANDAGES
Elastic retention bandage. It provides wound care splint retention and light support.
Confortable non woven fabric (cotton, viscose and polyamide) coated with syntethic latex.
Stretched length 4 m.

- BSN code: 45473-00
- Elastomull® Haft 4 m x 10 cm
- Box of 10 pcs.
- Multilanguage box: GB, FR, IT, ES, DE, PT, NL, SE.