



| | |
|---------------|----------------|
| Product Code | 59219 |
| Unit of sale | 1 pc |
| Minimum order | 1 |
| Type | Medical device |
| Class | I |

| | |
|------------|---------------|
| UK-REP | No |
| CH-REP | No |
| RDM (NSIS) | 2294821 |
| CND | M0303010201 |
| EAN/UPC | 4042809022063 |

Description

ELASTOMULL® COHESIVE BANDAGES

Elastic retention bandage. It provides wound care splint retention and light support.
Comfortable non woven fabric (cotton, viscose and polyamide) coated with syntethic latex.
Stretched length 20 m.

- BSN code: 45478-00
- Elastomull® Haft 20 m x 10 cm
- Box of 1 pc.
- Multilanguage box: GB, FR, IT, ES, DE, PT, NL, SE.