

PEHA-HAFT 300037 ELASTIC COHESIVE BANDAGE 4 m
x 2.5 cm - latex free



Product Code	59230
Unit of sale	box of 8 pcs
Minimum order	1
Type	Medical device
Class	I

RDM (NSIS)	2523552
CND	M0303010201
EAN/UPC	4049500030915
PARAF	924277837

Description

PEHA-HAFT® ELASTIC COHESIVE BANDAGES - LATEX FREE

- Hartmann code: 300037
- Peha-haft® 4 m x 2.5 cm
- Minimum order: box of 8

Single-use cohesive conforming bandages for retention of primary and secondary dressings (such as splints, casts, catheters and cannulas) particularly on joints as well as on conical and round parts of the body.
Double adhesive effect due to its crêpe texture and micro-structured coating.

- 85 % extensibility
- low consumption of material due to its strong self-adhesive effect
- bandage adheres to itself, but not to skin, hair or clothing
- permeable to air and skin-friendly.

Ideal for the stabilization of muscles, tendons, ligaments, joints with low immobilization effect. Also used in medical interventions when acute and chronic wounds need to be treated short or long term. Class I. Medical devices. It is made of 34% viscose, 34% cotton and 32% polyamide.

Multilingual box: GB, FR, IT, ES, PT, DE, NL, SE, FI, DK, PL, CZ, HU, SI, GR, Arabic.