

PEHA-HAFT 300044 ELASTIC COHESIVE BANDAGE 20
m x 6 cm - latex free



Product Code	59237
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

RDM (NSIS)	2523559
CND	M0303010201
EAN/UPC	4049500750042
PARAF	921568352

Description

PEHA-HAFT® ELASTIC COHESIVE BANDAGES - LATEX FREE

- Hartmann code: 300044
- Peha-haft® 20 m x 6 cm
- Minimum order: box of 1

Single-use cohesive conforming bandages for retention of primary and secondary dressings (such as splints, casts, catheters and cannulas) particularly on joints as well as on conical and round parts of the body.
Double adhesive effect due to its crêpe texture and micro-structured coating.

- 85 % extensibility
- low consumption of material due to its strong self-adhesive effect
- bandage adheres to itself, but not to skin, hair or clothing
- permeable to air and skin-friendly.

Ideal for the stabilization of muscles, tendons, ligaments, joints with low immobilization effect. Also used in medical interventions when acute and chronic wounds need to be treated short or long term. Class I. Medical devices. It is made of 34% viscose, 34% cotton and 32% polyamide.

Multilingual box: GB, FR, IT, ES, PT, DE, NL, SE, FI, DK, PL, CZ, HU, SI, GR, Arabic.