

CLOSED LAUNDRY TROLLEY - painted steel



Product Code	45925
Unit of sale	1 pc
Minimum order	1
Type	No medical device
Class	

UK-REP	No
CH-REP	No
EAN/UPC	8023279459258

Description

CLOSED LAUNDRY TROLLEY - painted steel

Trolleys for collection of dirty laundry and distribution of clean linen to the patient.
Stainless steel or painted steel frame.
Equipped with 4 Ø 125 mm anti-thread swiveling castors (2 with brake), 2 rectangular supports with rubber bag stop ring, 2 lids and 3 shelves of which 2 adjustable.

Suitable for 70 litres waste bags.

Size: 1330x650xh 1200 mm (45925) 1330x650xh 1140 mm (45926)

Made in Italy.