



Product Code	59300
Unit of sale	box of 250 pcs
Minimum order	1
Type	Medical device
Class	II A

UK-REP	Yes
--------	-----

RDM (NSIS)	2620360
CND	H900102
EAN/UPC	4049500263788
UMDNS	10-288

## Description

### OMNISTRIP® WOUND CLOSURE STRIPS - STERILE

- OMNISTRIP® - 3x76 mm
- Hartmann code: 540 681
- Minimum order: box of 250 (50 bags of 5)

Omnistrip® is a wound closure strip for atraumatic closure of minor and surgical incisions in primary healing wounds. It protects the wound margins against tension and stabilizes fresh scars, once sutures or staples have been removed. Rounded edges, conformable and stretchable, with a high porosity, permeability and breathability and good adherence to the skin. Radiolucent, temperature resistant and durable.

Made of a carded polyester-/polyamide skin-coloured non-woven, which is coated with an acrylic adhesive.

Multilingual manual and box: GB, FR, IT, ES, PT, DE, PL, NL, SE, DK, NO, FI, GR, CZ, SK, HU, RU, Arabic.