



Product Code	59301
Unit of sale	box of 300 pcs
Minimum order	1
Type	Medical device
Class	II A

UK-REP	Yes
RDM (NSIS)	2620361
CND	H900102
EAN/UPC	4049500263801
UMDNS	10-288
PARAF	987475492

## Description

### OMNISTRIP® WOUND CLOSURE STRIPS - STERILE

- OMNISTRIP® - 6x38 mm
- Hartmann code: 540 682
- Minimum order: box of 300 (50 bags of 6)

Omnistrip® is a wound closure strip for atraumatic closure of minor and surgical incisions in primary healing wounds. It protects the wound margins against tension and stabilizes fresh scars, once sutures or staples have been removed. Rounded edges, conformable and stretchable, with a high porosity, permeability and breathability and good adherence to the skin. Radiolucent, temperature resistant and durable.

Made of a carded polyester-/polyamide skin-coloured non-woven, which is coated with an acrylic adhesive.

Multilingual manual and box: GB, FR, IT, ES, PT, DE, PL, NL, SE, DK, NO, FI, GR, CZ, SK, HU, RU, Arabic.