



Product Code	59303
Unit of sale	box of 500 pcs
Minimum order	1
Type	Medical device
Class	II A

UK-REP	Yes
RDM (NSIS)	2620363
CND	H900102
EAN/UPC	4049500263849
UMDNS	10-288
PARAF	987475480

Description

OMNISTRIP® WOUND CLOSURE STRIPS - STERILE

- OMNISTRIP® - 6x101 mm
- Hartmann code: 540 684
- Minimum order: box of 500 (50 bags of 2x5)

Omnistrip® is a wound closure strip for atraumatic closure of minor and surgical incisions in primary healing wounds. It protects the wound margins against tension and stabilizes fresh scars, once sutures or staples have been removed. Rounded edges, conformable and stretchable, with a high porosity, permeability and breathability and good adherence to the skin. Radiolucent, temperature resistant and durable.

Made of a carded polyester-/polyamide skin-coloured non-woven, which is coated with an acrylic adhesive.

Multilingual manual and box: GB, FR, IT, ES, PT, DE, PL, NL, SE, DK, NO, FI, GR, CZ, SK, HU, RU, Arabic.