



|               |                   |
|---------------|-------------------|
| Product Code  | 35936             |
| Unit of sale  | 1 pc              |
| Minimum order | 1                 |
| Type          | No medical device |
| Class         |                   |

|         |               |
|---------|---------------|
| UK-REP  | Yes           |
| CH-REP  | Yes           |
| EAN/UPC | 8023279359367 |

## Description

### EURONDA EUROSEAL® INFINITY SEALING MACHINE

Euroseal® Infinity is a practical, sturdy, reliable and easy to use electronic thermosealing machine for continuous cycle processing at constant temperature.

#### Compact design

Compact and functional design to minimise the space required for its use so it can be comfortably placed on a table top or wall-mounted (need optional kit) in a vertical position.

#### Versatile

The roll holder is adjustable in two positions: high for best performance when the machine is wall-mounted, or low to speed up roll selection and changing operations.

#### Ergonomic

The new built-in cutter and sealing lever positioned at the front, not only improve ergonomics and safety, but also make packaging and sealing operations easier.

#### Safe

It has a metal bar that holds and locks the roll and ensures that it remains inside the sealing area. The smooth surfaces make it easy to clean.

#### Reliability

The quality of the materials and sturdy design make it indestructible, guaranteeing long-lasting efficiency and reliability.

Compliant with UNI EN 868-5 and EN 11607-2.

Multilingual manual: GB, IT.

Available QR code manual in: FR, ES, PT, DE, PL, CZ, SI, RO, RS, RU.

## Technical Specifications

Sealing size: 12x310 mm

Sealing temperature: 170-200°C

Size - weight: 430x430xh 215 mm - 7.1 kg

Power supply: 200-240 V; 50/60 Hz, 100W