

LAPARATOMIC SPONGE GAUZE 20x20 cm 4 plies



Product Code	59391
Unit of sale	box of 600 pcs
Minimum order	1
Type	Medical device
Class	II A

RDM (NSIS)	50151
CND	M0201030202
EAN/UPC	8030413011449
PARAF	987472204

Description

LAPAROTOMIC SPONGE GAUZE - NOT STERILE

- Laparotomic sponge gauze size: 20x20 cm
- Plies: 4
- Threads: 20

100% cotton gauze made of 4 layers (8 layers for 59394) of 20 threads/cm², sewn on sides and highly absorbent. Barium X-ray thread woven into the gauze, made of biocompatible polypropylene, with 60% barium sulphate. Class IIa Medical Device according to 93/42/CEE, 2007/47/CEE
Latex free.