

LAPARATOMIC SPONGE GAUZE 50x50 cm 4 plies



Product Code	59395
Unit of sale	box of 250 pcs
Minimum order	1
Type	Medical device
Class	II A

RDM (NSIS)	50151
CND	M0201030202
EAN/UPC	8030413003192
PARAF	987475124

Description

LAPAROTOMIC SPONGE GAUZE - NOT STERILE

- Laparotomic sponge gauze size: 50x50 cm
- Plies: 4
- Threads: 20

100% cotton gauze made of 4 layers (8 layers for 59394) of 20 threads/cm², sewn on sides and highly absorbent. Barium X-ray thread woven into the gauze, made of biocompatible polypropylene, with 60% barium sulphate. Class IIa Medical Device according to 93/42/CEE, 2007/47/CEE. Latex free.