

COTTON GAUZE ROLLS 80 cm x 1 m



Product Code	59410
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	Yes
--------	-----

EAN/UPC	5605622400143
---------	---------------

Description

100% COTTON GAUZE ROLLS

100% cotton gauze rolls with 20 threads sq cm according to European Pharmacopoeia and EN 14079.
X-Ray detectable thread (approx. 60% barium sulphat) 80 cm width, folded 8 to 10 cm ply.

- Extended width: 80 cm x 1 m
- Folded width: 10 cm x 1 m
- Plies: 8
- Threads sq cm: 20

Multilanguage box: GB, FR, IT, ES, PT, DE, NL, SE.