

COTTON RIBBON GAUZE 5 m x 1.25 cm - sterile



Product Code	59440
Unit of sale	box of 60 pcs
Minimum order	1
Type	Medical device
Class	II A

UK-REP	Yes
--------	-----

EAN/UPC	8023279594409
---------	---------------

Description

100% COTTON RIBBON GAUZE weight 33 g/sq. m

22 threads/sq. m

5 m x 1.25 cm - sterile

- Latex and PVC free
- Absorption of fluids inside cavities
- Wound care during surgery

Multilanguage box: GB, FR, IT, ES, PT, DE, NL, SE.