

COTTON RIBBON GAUZE 5 m x 2.5 cm - sterile



Product Code	59441
Unit of sale	box of 30 pcs
Minimum order	1
Type	Medical device
Class	II A

UK-REP	Yes
--------	-----

EAN/UPC	8023279594416
---------	---------------

Description

100% COTTON RIBBON GAUZE weight 33 g/sq. m

22 threads/sq. m

5 m x 2.5 cm - sterile

- Latex and PVC free
- Absorption of fluids inside cavities
- Wound care during surgery

Multilanguage box: GB, FR, IT, ES, PT, DE, NL, SE.