


BOWIE DICK TEST PACK



Product Code	35950
Unit of sale	box of 20 pcs
Minimum order	1
Type	No medical device
Class	

UK-REP	Not required
CH-REP	Not required

EAN/UPC	8023279359503
---------	---------------

Description

BOWIE DICK TEST PACK

- box of 20 pcs
- dimensions: 123 x 123 x 16 mm

Used for daily monitoring of pre-vacuum steam sterilizer.
Based upon the Bowie Dick test result the steam sterilizer is released for usage during the day.
The Bowie Dick test detect failures like air leaks, insufficient vacuum, inadequate steam penetration.
Package and manual in: GB, IT, ES, FR, PT, GR, DE.

- fulfil II requirements of EN 285 and ISO 17665
- conform to ISO 11140 - 4 Type 2 / EN 867 - 4 class B (EU standard 7 kilo)
- independently tested to show performance similar to original pack
- clear blue - violet colour change
- easy to interpret and document
- non-Toxic / lead free