


BOWIE DICK HELIX TEST PACK



Product Code35951

Unit of salebox of 100 pcs

Minimum order1

TypeNo medical device

Class

UK-REPNot required

CH-REPNot required

EAN/UPC8023279359510

Description

BOWIE DICK HELIX TEST PACK

- box of 100 pcs
- includes: cotton bag, helix device (tube 1.5 m), 100 tests and user manual GB, FR, IT, ES

Helix devices and their indicators have proven that they detect failures during each cycle much earlier than other chemical indicators / integrators.

The EN 285: 2009 are specifying the Helix test (Hollow A) to be used as Bowie Dick Test as well as Batch Monitoring when hospitals are sterilizing hollow loads, next to that this Bowie Dick Helix Test is an obligatory test for class B Bench Top Steam sterilizers.

- conform to ISO 11140 - 1 Type 2 - EN 867 part 5 class B
- clear transition of color into pink
- easy to interpret and document
- non-Toxic / Lead Free
- PCD guaranteed for 500 Cycles (134°C - 3,5 min)