

HYPAFIX DRESSING RETENTION 10 m x 50 mm



Product Code	59510
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

RDM (NSIS)	2265105
CND	M050201
EAN/UPC	4042809000702

Description

BSN HYPAFIX®

Non woven dressing retention sheet for the retention of dressings especially on highly contoured and frequently moved parts of the body. Conformable non woven fabric coated with a poliacyrylic adhesive mass.

HYPAFIX® 50 mm x 10 m
BSN code: 71443-01