

OMNIFIX ELASTIC HARTMANN NON WOVEN RETENTION TAPE 10 m x 10 cm - 900603



Product Code	59522
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

RDM (NSIS)	2454298
EAN/UPC	4049500967778
UMDNS	10-028
GMDN	16866

Description

OMNIFIX® ELASTIC NON-WOVEN RETENTION TAPE - STERILE
- OMNIFIX® elastic - 10 m x 10 cm
- HARTMANN code: 900603

Omnifix® Elastic is a hypoallergenic non-woven fibre fleece fixation tape with synthetic adhesive that offers comfortable, full-coverage bandage retention wrap that can be stretched laterally. It is soft and elastic as well as breathable and permeable to water vapor reducing the risk of maceration.

-skin friendly, synthetic adhesive that is suitable for patients with delicate skin
-easy to apply, a measuring scale printed on the cover paper facilitates exact trimming; the wave-shaped, die-cut backing paper is easy to remove
-easy to remove painlessly and without leaving residues
-secure adhesion

Applied as a dressing retention by full covering, especially on joints and angular parts of the body without constraining movement; for securing equipment, probes and cannulae.

Multilingual box: GB, FR, IT, ES, PT, DE, NL, DK, NO, PL, CZ, SK, HU, SI, HR, RO, BG, GR, RU, TR, CN, Arabic.