

# OMNIFIX ELASTIC HARTMANN NON WOVEN RETENTION TAPE 10 m x 20 cm - 900605



Product Code	59524
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

RDM (NSIS)	2454300
EAN/UPC	4049500967792
UMDNS	10-028
GMDN	16866

## Description

OMNIFIX® ELASTIC NON-WOVEN RETENTION TAPE - STERILE  
- OMNIFIX® elastic - 10 m x 20 cm  
HARTMANN code: 900605

Omnifix® Elastic is a hypoallergenic non-woven fibre fleece fixation tape with synthetic adhesive that offers comfortable, full-coverage bandage retention wrap that can be stretched laterally. It is soft and elastic as well as breathable and permeable to water vapor reducing the risk of maceration.

skin friendly, synthetic adhesive that is suitable for patients with delicate skin  
easy to apply, a measuring scale printed on the cover paper facilitates exact trimming; the wave-shaped, die-cut backing paper is easy to remove  
easy to remove painlessly and without leaving residues  
secure adhesion

Applied as a dressing retention by full covering, especially on joints and angular parts of the body without constraining movement; for securing equipment, probes and cannulae.

Multilingual box: GB, FR, IT, ES, PT, DE, NL, DK, NO, PL, CZ, SK, HU, SI, HR, RO, BG, GR, RU, TR, CN, Arabic.