

OMNIFIX ELASTIC HARTMANN NON WOVEN RETENTION TAPE 10 m x 30 cm - 900607



| | |
|---------------|----------------|
| Product Code | 59525 |
| Unit of sale | 1 pc |
| Minimum order | 1 |
| Type | Medical device |
| Class | I |

| | |
|------------|---------------|
| RDM (NSIS) | 2454301 |
| EAN/UPC | 4049500203326 |
| UMDNS | 10-028 |
| GMDN | 16866 |

Description

OMNIFIX® ELASTIC NON-WOVEN RETENTION TAPE - STERILE - OMNIFIX® elastic - 10 m x 30 cm - HARTMANN code: 900607

Omnifix® Elastic is a hypoallergenic non-woven fibre fleece fixation tape with synthetic adhesive that offers comfortable, full-coverage bandage retention wrap that can be stretched laterally. It is soft and elastic as well as breathable and permeable to water vapor reducing the risk of maceration.

skin friendly, synthetic adhesive that is suitable for patients with delicate skin - easy to apply, a measuring scale printed on the cover paper facilitates exact trimming; the wave-shaped, die-cut backing paper is easy to remove - easy to remove painlessly and without leaving residues - secure adhesion

Applied as a dressing retention by full covering, especially on joints and angular parts of the body without constraining movement; for securing equipment, probes and cannulae.

Multilingual box: GB, FR, IT, ES, PT, DE, NL, DK, NO, PL, CZ, SK, HU, SI, HR, RO, BG, GR, RU, TR, CN, Arabic.