

SILVERCEL 3M NON-ADHERENT HYDRO-ALGINATE ANTIMICROBIAL DRESSING WITH SILVER 11x11 cm



Product Code	59571
Unit of sale	box of 10 pcs
Minimum order	1
Type	Medical device
Class	III
RDM (NSIS)	2376702
CND	M04040802

Description

3M™ SILVERCEL™

3M CODE: CAD7011

3M™ Silvercel™ Non-Adherent Antimicrobial Dressing with Silver 11x11 cm - box of 10

3M™ Silvercel™ dressing is a sterile, non-woven pad composed of alginate, carboxymethylcellulose (CMC), and silver-coated nylon fibers that help maintain a moist wound environment. The excellent fluid absorption ability and lasting antimicrobial activity make this dressing ideal for chronic wounds with moderate to high levels of exudate, such as pressure ulcers, venous ulcers, neuropathic (diabetic) ulcers, donor sites, post-operative wounds, and trauma wounds.

The antimicrobial effect is secured by a sustained release of silver ions, which can last up to 7 days. This dressing acts against a broad spectrum of germs that are pathogenic to the wound, including MRSA, MRSE, and VRE.

Easylift™ Technology

Easylift™ Precision Film Technology is built into the non-adherent wound contact layer. It enables pain-free removal and provides effective protection of newly-formed tissue while minimizing fiber detachment.