



Product Code	59575
Unit of sale	box of 25 pcs
Minimum order	1
Type	Medical device
Class	III

RDM (NSIS)	2487945
CND	M040412
PARAF	901685279
GTIN	04054596954595
GMDN	47203

## Description

### INADINE™ PVP-I NON ADHESIVE DRESSING

- Inadine™ 5x5 cm

Non adherent dressing, made of low adherent knitted viscose fabric impregnated with a polyethylene glycol (PEG) base, containing 10% Povidone Iodine.

Designed to protect wounds even when they are infected, providing long-lasting antimicrobial effect against bacteria, mycobacteria, fungi, protozoa and viruses, including resistant species such as MRSA.

Inadine dressing minimizes adherence to the wound bed, therefore reducing pain for patient and the risk of damage to the granulation tissue at dressing removal.

Colour evolution, from orange to white indicates loss of antiseptic efficacy and dressing needs to be changed.

Weight: 3.25g.

Indicated for ulcerative wounds, prevention of infection in minor burns and traumatic skin loss injures and heavily infected wounds.

Latex free.

Individually packed in aluminium pouch.

Gamma ray sterilization.