



Product Code	59576
Unit of sale	box of 10 pcs
Minimum order	1
Type	Medical device
Class	III

RDM (NSIS)	2487946
CND	M040412
PARAF	901993651
GTIN	04054596954625
GMDN	47203

Description

INADINE™ PVP-I NON ADHESIVE DRESSING

- Inadine™ 9.5x9.5 cm

Non adherent dressing, made of low adherent knitted viscose fabric impregnated with a polyethylene glycol (PEG) base, containing 10% Povidone Iodine.

Designed to protect wounds even when they are infected, providing long-lasting antimicrobial effect against bacteria, mycobacteria, fungi, protozoa and viruses, including resistant species such as MRSA.

Inadine dressing minimizes adherence to the wound bed, therefore reducing pain for patient and the risk of damage to the granulation tissue at dressing removal.

Colour evolution, from orange to white indicates loss of antiseptic efficacy and dressing needs to be changed.

Weight: 3.25g.

Indicated for ulcerative wounds, prevention of infection in minor burns and traumatic skin loss injures and heavily infected wounds.

Latex free.

Individually packed in aluminium pouch.

Gamma ray sterilization.