



| | |
|---------------|----------------|
| Product Code | 59580 |
| Unit of sale | box of 10 pcs |
| Minimum order | 1 |
| Type | Medical device |
| Class | II B |

| | |
|------------|---------------|
| RDM (NSIS) | 2477443 |
| CND | M040299 |
| EAN/UPC | 4054596952850 |
| PARAF | 983199062 |

Description

KERRAMAX CARE™ SUPER-ABSORBENT DRESSING

3M™ code: PRD500-050
 Kerramax™ Care 10x10 cm
 Minimum order: Box of 10

The unique wicking layer distributes fluid both horizontally and vertically throughout the dressing, utilizing the full absorption capacity. Exu-Safe™ Technology, built into the dressing core, is an advanced method of absorbing high level of exudate from wounds, including bacteria and MMPs, creating an optimal environment for wound healing. It protects delicate wound tissue and surrounding, helping to reduce the risk of maceration, even under compression. Its heat-sealed border prevents exudate leakage.

It can be left in place for up to 7 days.
 Absorption capacity: 0.95 g/cm².
 Indicated for: pressure injuries venous leg ulcers and diabetic foot ulcers.

Latex free.
 Ethylene oxide sterilization.